

Mobile device for managing medical procedures with Closed technical compartment with a high degree of sanitisation.



SA.F.E.




- Suitable for Operating Department, Intensive Care Unit, Inpatient Department.
- Frames for ISO baskets*
- Made from Austenitic Stainless Steel AISI 304
- Work Surface:
640 x 460 mm. (Larg. x Prof.)
- No. 2 directional handles.
- No. 4 Twin Conductive Wheels Ø 125 mm. (No. 2 with Brake)
- Perimeter Bumpers
- Maximum capacity 2,5 STU / 2,5 ISO
- Overall dimensions
821 x 621 x 1015 mm.
(W x D x H)
SA.F.E. Tech.


*Baskets NOT INCLUDED but optional:

- ABS
- AISI 304
- AISI 304 Wire

Product Highlights




Sanitizable Handle




*Edges/Corners internal
R=20 mm*



Exhaust manifold Macroimpurity



*Colors/treatments
For identification and detection of
macroimpurity*



**Patented
Device**



Great Handling



**Washing also
with thermodisinfectant**

Product Highlights (Frames for ISO baskets)



ADVANTAGES

- Interchange of baskets between-FRENDY-trolleys-Modular Systems.
- Greater Sanifiability (total extraction for sanitization or sterilization)



Format: A1

- N°3 FRAMES (removables) in Austenitic Stainless Steel (AISI 304).
(W x D x H) 600 x 400 x 150 mm.
- N°1 FRAME (removables) in Austenitic Stainless Steel (AISI 304).
(W x D x H) 600 x 400 x 300 mm.



Format: **A2**

- N°5 FRAMES (removables) in Austenitic Stainless Steel (AISI 304).
(W x D x H) 600 x 400 x 150 mm.



Format: **A3**

- N°6 FRAMES (removables) in Austenitic Stainless Steel (AISI 304).
(W x D x H) 300 x 400 x 150 mm.
- N°1 FRAME (removables) in Austenitic Stainless Steel (AISI 304).
(W x D x H) 600 x 400 x 300 mm.



Format: **A4**


- ABS Panels for ISO Baskets 600 x 400
- Bifront access


(Emergency)





Format: **E1 (Emergency)**






- N°3 FRAMES (removables) in Austenitic Stainless Steel (AISI 304).
(W x D x H) 600 x 400 x 150 mm.
- N°1 FRAMES (removables) in Austenitic Stainless Steel (AISI 304).
(W x D x H) 600 x 400 x 300 mm.
- N°3 ISO Baskets (ABS) H=100
- N°1 ISO Basket (ABS) H=200
- I.V. Stand
- Table for CPR & supports
- Aspirator Support
- Oxygen Tank Holder
- Swivel support for defibrillator
- Opening device with security seal (Tamper-proof divider)





FR-ACS00	Removable Side Shelf Capacity A4. Usable support dimension (W X D) 320 x 320 mm (+8/-0)		172
FR-ACS01	Complete accessories support		173
FR-ACS02	Partial accessories support		174
COS 6150	Dispenser for various materials. Made entirely from aAustenitic Stainless Steel AISI 304. No. 5 + 4 non-toxic plastic containers. Overall Dimensions (W x D x H) 620 x 180 x 420 mm		311





FR-ACS03	Swivel support for defibrillator		175
FR-ACS04	Aspirator Support		176
FR-ACS05	Oxygen Tank Holder		177
FR-ACS06	Table for CPR & supports		178
FR-ACS07	Opening device with security seal (Tamper-proof divider)		179

FR-ACS08	Key central opening and closing device		180
FR-ACS09	Waste bin with knee opening		181
FR-ACS10	Waste bin for sharp/dangerous content		182
FR-ACS11	Perforated plate probe basket 120 x 170 x 600H (ACS 570-1720 with lateral restraint)		183
FR-ACS12	Electrification through no. 4 electrical sockets and bipolar switch (CE Conformity 93/42)		184

FR-ACS 13	External Bottle Holder (No. 4 for 500 ml)		185
FR-ACS14	External Holder for A4 Gatherers (No. 3 High 35 mm)		186

ISO 6410X	ISO basket. Made from Austenitic Stainless Steel AISI 304. Overall dimensions (W x D x H) 600x400x100 mm		154
ISO 6420X	Basket ISO. Made entirely from Austenitic Stainless Steel AISI 304 Overall dimensions (W x D x H) 600x400x200 mm		155
ISO 6405S	ISO basket. Made from ABS. Overall dimensions (W x D x H) 600x400x50 mm		158
OGDZ20204A	Partition (2X2) 4 Compartments for Basket ISO 6405S		159
OGDZ20306A	Partition (2X3) 6 Compartments for Basket ISO 6405S		160
OGDZ30412A	Partition (3X4) 12 Compartments for Basket ISO 6405 S		161

ISO 6410S	ISO basket. Made from ABS. Overall dimensions (W x D x H) 600x400x100 mm		162
OGDZ20204B	Partition (2X2) 4 Compartments for Basket ISO 6410S		163
OGDZ20306B	Partition (2X3) 6 Compartments for Basket ISO 6410S		164
OGDZ30412B	Partition (3X4) 12 Compartments for Basket ISO 6410S		165
ISO 6420S	ISO basket. Made from ABS. Overall dimensions (W x D x H) 600x400x200 mm		166
OGDZ20204C	Partition (2X2) 4 Compartments for Basket ISO 6420S		167

OGDZ20306C	Partition (2X3) 6 Compartments for Basket ISO 6420S		146
OGDZ30412C	Partition (3X4) 12 Compartments for Basket ISO 6420S		147
ISO 3410X	Basket ISO. Made entirely from Austenitic Stainless Steel AISI 304 Overall dimensions (W x D x H) 300x400x100 mm		156
ISO 3420X	Basket ISO. Made entirely from Austenitic Stainless Steel AISI 304 Overall dimensions (W x D x H) 300x400x200 mm		157
ISO 6420TX	ISO Basket in steel rod. Made entirely from Austenitic Stainless Steel AISI 304. Overall dimensions (W x D x H) 600x400x180 mm		170



**PROGETTAZIONE, PRODUZIONE, VENDITA E ASSISTENZA RELATIVA A:
DISPOSITIVI PER REPARTI OPERATORI ED AMBIENTI AD ELEVATO GRADO DI SANIFICABILITÀ.
DISPOSITIVI PER LAPAROSCOPIA**

- Trasporto, Stoccaggio e Gestione di Dispositivi medici e non.



- Dispositivi Riscaldatori termoregolati ad Aria o Acqua.



- Dispositivi per la Gestione dei Dispositivi Medici e Presidi



- Dispositivi In assistenza agli Operatori



- Dispositivi Elevatori di Soluzione Infusionale



- Sedili per Operatori



- Dispositivi Integrativi dell'ergonomia (Pedane, Tutori per Operatori)



- Vassoi e Bacinelle ad uso Clinico.



- Supporti Unificati per la gestione dei Dispositivi Medici e attrezzature nei Reparti



- Postazioni per sanificazione Operatori di Camera Operatoria (Lavelli-Lavamani ed accessori)



- LaparoTensor®



LUCINI SURGICAL CONCEPT srl
Via Sant'Alberto, 4/6/8
20020 Lainate - Milano - Italy
Tel. +39 (0)2 91 86 543 - Fax +39 (0)2 91 89 185
E-mail: info@lucinisurgicalconcept.com
<http://www.lucinisurgicalconcept.com>

